



“AN
ISO 9001- 2015
WHO - GMP
CERTIFIED CO.”



- * HOSPITAL FURNITURE
- * MEDICAL EQUIPMENTS
- * SCIENTIFIC INSTRUMENTS
- * HOSPITAL HOLLOWWARE

ICU Bed Five Functional Electric



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest, trendelenburg/Reverse trendelenburg hi-low position obtained by electric actuator
- SS molded head and foot panels
- Aluminium collapsible side safety railing
- Bed mounted on 5 inch special casters, two with brakes
- One S.S.I.V rod
- Overall size: 220 L * 100 W * 60-80 H cm
- Pre-treated and Epoxy powder coated



ICU Bed Five Functional Electric (Eco Model)



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest, trendelenburg/Reverse trendelenburg hi-low position obtained by electric actuator
- ABS molded head and foot panels
- Aluminium collapsible side safety railing
- Bed mounted on 5 inch special casters, two with brakes
- One S.S.I.V rod
- Overall size: 220 L * 100 W * 60-80 H cm
- Pre-treated and Epoxy powder coated

ICU Bed Five Functional Manual



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest, trendelenburg/Reverse trendelenburg hi-low position obtained by crank mechanism
- ABS molded head and foot panels
- ABS side safety railing
- Bed mounted on 5 inch special casters, two with brakes
- One S.S.I.V rod
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated



ICU Bed Five Functional Manual (Eco Model)



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest, trendelenburg/Reverse trendelenburg hi-low position obtained by crank mechanism
- ABS molded head and foot panels
- Aluminium collapsible side safety railing
- Bed mounted on 5 inch special casters, two with brakes
- One S.S.I.V rod
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest operated by Motor & Remote
- ABS molded head and foot panels
- Aluminium collapsible side safety railing
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated
- Bed mounted on 5 inch special casters, two with brakes

Fowler Bed Electrical



Fowler Bed Super Deluxe



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest obtained by crank mechanism
- ABS molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated
- Bed mounted on 5 inch special casters, two with brakes



Fowler Bed Deluxe A



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest obtained by crank mechanism
- Full SS molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated

- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest obtained by crank mechanism
- SS cum sunmica molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated

Fowler Bed Eco



Fowler Bed Deluxe



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest, knee-rest obtained by crank mechanism
- Powder Coated molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated



Semi Fowler Bed Electric



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest obtained by Motor & Remote
- ABS molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated

Semi Fowler Bed Super Deluxe



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest obtained by crank mechanism
- ABS molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated

Semi Fowler Bed Deluxe



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest obtained by crank mechanism
- SS cum Sunmica molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated



Semi Fowler Bed Eco Model



- Frame made of strong rectangular CRCA pipe
- Uniformly perforated Four Sectional CRCA Sheet top
- Back-rest obtained by crank mechanism
- Powder Coated molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated

Plain Bed Super Deluxe



- Frame made of strong rectangular CRCA pipe
- ABS molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated

Plain Bed Deluxe



- Frame made of strong rectangular CRCA pipe
- SS molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated



Plain Bed Eco Model



- Frame made of strong rectangular CRCA pipe
- Powder Coated molded head and foot panels
- Overall size: 220 L * 100 W * 60-80 H cm
- Pretreated and Epoxy powder coated

Pediatric Bed



- Main frame of rectangular pipe
- Perforated CRCA sheet top.
- Tubular head & foot bows of equal height with vertical tube supports on PVC stumps.
- Drop down type safety side railings are provided.
- Size: 137L x 76W x 60H cms.
- Finish Pre-treated & Epoxy Powder Coating finish.

Baby Bassinet



- Measurement: 870x530x780~980mm
- Aluminum column, height is adjustable
- Weighing system is optional



O.T Table Electric



- SS five sectional top with large cut for perianal drainage tray, built Kidney Bridge
- Both side Controls
- Electrically Operated
- Heavy duty steel base with nylon wheels for easy mobility or stability with single action foot actuated floor locking
- Smooth gear mechanism to avail all operating positions
- 50, Thick electrically sealed mattress
- Provide with all sporting accessories
- Durable heavy duty leak proof Motors & Remote
- Provided with standard accessories

OT Table Hydraulic




- SS five sectional top with large cut for perianal drainage tray, built Kidney Bridge
- Both side Controls
- Hydraulic Operation tables
- Heavy duty steel base with nylon wheels for easy mobility or stability with single action foot actuated floor locking
- Smooth gear mechanism to avail all operating positions
- 50, Thick electrically sealed mattress
- Provide with all sporting accessories
- Durable heavy duty leak proof hydraulic pump with double cylinder concealed in central chamber
- Provided with standard accessories

O.T Table C-Arm Electric



- Five Sectional Stainless Steel Radio Translucent Top.
- Eccentrically Positioned top for maximum access for movement of c-arm image intensifier.
- Base and cover made of heavy & hygienic stainless steel sheet.
- Mobile on Smooth running castors and floor locking system for stability.
- Remote Control system Offers Define Position:
 - Up and Down Position
 - Chair position for Seating operation

O.T Table C-Arm Hydraulic



- Five Sectional Stainless Steel Radio Translucent Top.
- Eccentrically Positioned top for maximum access for movement of c-arm image intensifier.
- Base and cover made of heavy & hygienic stainless steel sheet.
- Mobile on Smooth running castors and floor locking system for stability.
- Durable heavy duty leak proof hydraulic pump with double cylinder concealed in central chamber

Delivery Table



- Overall approx. size: 187 L * 51 W * 81 H
- Adjustable backrest by gas spring system
- Two section foam padded, upholstered top
- Three drawers
- Three cabinets
- Tray for BP apparatus
- Pretreated and Epoxy powder coated
- Adjustable lithotomy poles with straps
- Also available in SS

Delivery Table Three Fold



- Made In Three Sections, Adjustable back rest & Leg Section With lever mechanism to various Positions capable of tilting for obstetric and gynecological Use
- Epoxy Coated Tubular Steel frame
- Highly Coated mattress pad
- Supplied with two leg holders with canvas strap

Delivery Bed



- Tubular Frame work mounted on 125mm dia castor two with brakes
- Three Section top made of MS Pipe with fix cushion mattress
- Height adjustment & back section adjustment through electric actuator system
- Leg end section can slide under the main section
- SS telescopic IV Rod

Delivery Table Hydraulic Deluxe

- Three Sectional Stainless Steel Top, Back-rest adjustable by Ratchet
- Middle Section with "U" cut, leg section can slide under main frame
- Trendelenburg Position achieved with Pneumatic pump system
- Provide with head & side safety rails, SS Hand Grips
- Three sectional Matters, a pair of Lithotomy crutches & one IV Rod
- Sliding SS Douche Tray
- Legs Fitted on four rubber stumps
- Overall Extendable Size: 187 L x 96 W x 80 H C



Delivery Table Electric



- The Table is design for women's obstetric & Gynecological Operation
- Operation Position - up / Down, adjustable back rest and reverse Position
- Main frame work Epoxy Powder coated Steel
- Supplied with two leg holders with canvas strap
- Adjustable legs support & Sliding Grip Holder
- Equipped with 8cm Thick High Quality Mattress

Mobile O.T Light

- Have daylight characters
- can be adjustable as per the operation requirements
- lights arms are made of stainless steel fitting
- Range of the color temp $4200\text{ K} \pm 300$ as natural color tone
- Detachable handle for simplistic sterilization
- Diachronic reflectors & color correcting filters
- Metallic reflectors available



O.T Light Single Dome Halogen



- lux intensity 90,000
- halogen bulb of 24v 150w
- 20" dome diameter
- Arms And all Accessories are made up of SS

O.T Twin Light Halogen

- lux intensity 1,80,000
- halogen bulb of 24v 150w
- 20" dome diameter each
- Arms And all Accessories are made up of SS
- Also available in diachloride & Aluminum Domes



LED O.T Light Single Dome



- **intensity:120000 lux:10%**
- **size of led field:12-30 cm**
- **Number of Led:19x4 (multiple colour)**
- **With Brightness Controller (Feather Touch)**
- **Adjustable Focus**

LED O.T Light Single Dome



- **intensity:135000 +- 10%**
- **Number of LEDs :-51**
- **Size of led field : 18-35cm**
- **Adjustable Focus**
- **With anti hms system**
- **With feather touch controller**

LED O.T Light Twin Light



- **intensity:2,20,000 lux**
- **size of led field:12-30 cm**
- **Number of Led:19x8 (multiple colour)**
- **With Brightness Controller (Feather Touch)**
- **Adjustable Focus**

Stretcher Trolley



Stretcher Trolley Deluxe



Stretcher Trolley Cum Chair



Emergency Patient Trolley



Ambulance Stretcher



Scrub Station With Sensor / Without Sensor



Attendant Bed Cum Chair



Examination Table Plain



Blood Donating Chair



Hospital Three Fold Screen



Crash Cart With ABS Drawers



Crash Cart With S.S Drawers



Three Seater Chrome



Three Seater



Examination Couch Deluxe



Cylinder Trolley



Revolving Stool Super Deluxe



Revolving Stool Deluxe



Bed Side Locker S.S



Bed Side Locker Semi Deluxe



Bed Side Locker Deluxe



Bed Side Locker ABS Type



ABS Anesthetic Cart



Wheel Chair



Instrument Trolley All SS



Instrument Trolley Box Type



Instrument Trolley Powder Coated



Instrument Trolley Three Shelf



Dressing Trolley All S.S



E.C.G Trolley



Linen Trolley



Laparoscopy Trolley



Food Table



Bio-Waste Trolley



Mayo's Trolley



Instrument Cabinet




Spine Board



Scoop Stretcher



Tow Fold Stretcher



Four Fold Stretcher




I.V Stand



Kick Bucket



Oxygen Pipeline



Fogger Machine




Three Para Monitor



Five Para Monitor



E.C.G Machine



**New & Pre-Owned
Defibrillator**



**Cautery Machine
(300w / 400w / Sealer)**



Syringe Pump



**Fetal Monitors &
C.T.G Machines**



**Bi-pap & C-pap
Machines**



Oxygen Concentrators



**New & Pre-Owned
C.T Scan Machines**



**New & Pre-Owned
M.R.I Scan Machines**



**New & Pre-Owned
C.T Scan Machines**



**New & Pre-Owned
C.R Systems**



**New & Pre-Owned
C-arm Machines**



**New & Pre-Owned
Diagnostic Ultrasound Scanners**



**T.C & General
Instruments**



**Laparoscopy & Endoscopy
Equipments & Instruments**



**New & Pre-Owned
Endoscopy & Laparoscopy
Camera Units**



**New & Pre-Owned
Gastroscopy Units**



**New & Pre-Owned
Bubble C-pap**



Dental Chair



**Microscopes, Centrifuge,
Incubator, Hot Air Oven
& All Scientific Goods**



Hematology Analyzer



Electrolyte Analyzer



**Arterial Blood Gas
Machine**



Biochemistry Analyzer




Elisa Reader & Washer



Baby Radiant Warmer



LED Photo Therapy Machine
Single Surface / Double Surface



Baby Incubator



Horizontal & Vertical
Autoclave



E.T.O Sterilizer



Suction Machines



Anesthesia Workstation
& Trolleys



Ventilator Machine



Automatic External
Defibrillator



Bed Headed Panels



O.T Pendant



**E.N.T & Neuro
Microscopes**



Physiotherapy Goods



General Surgical Goods



**I.C.U & Wards
Curtaining**



**Ortho & Neuro
Attachments**



**Mannequins & All
Dummies**



**Dressing Drum, Kidney Tray
& Instrument Tray**





Surgical Tray with cover, Stainless Steel Seamless

Premium		Economy	
AHW-108	Size 200x150x50 mm (8x6x2 in)	AHWE-104	Size (9x6x2 in)
AHW-109	Size 250x200x50 mm (10x8x2 in)	AHWE-105	Size (8x6x2 in)
AHW-110	Size 300x250x50 mm (12x10x2 in)	AHWE-106	Size (10x8x2 in)
AHW-111	Size 350x250x50 mm (14x10x2 in)	AHWE-107	Size (11x7x2 in)
AHW-155	Size 310x195x63 mm (12.5x7.75x2.5 in)	AHWE-108	Size (12x8x2 in)
AHW-156	Size 225x125x50 mm (9x5x2 in)	AHWE-109	Size (12x10x2 in)
		AHWE-110	Size (14x10x2 in)
		AHWE-111	Size (15x12x2.5 in)
AHW-162	Surgical Tray without cover, Stainless Steel Seamless (Mayo Tray) Size 480x330x19 mm (19x13x0.75 in)	AHWE-112	Size (18x12x3 in)

Kidney Tray Stainless Steel

Premium		Economy	
AHW-113	Size (150x70x40 mm) 6"	AHWE-113	6"
AHW-114	Size (200x90x45 mm) 8"	AHWE-114	8"
AHW-115	Size (250x100x52 mm) 10"	AHWE-115	10"
AHW-116	Size (300x130x60 mm) 12"	AHWE-116	12"
AHW-157	250 X 150 X 40 (Cap. 825 ml) 10"		



Catheter Tray Stainless Steel

AHW-117	222x82x41 mm (8.5x3.25x1.5 in)
AHWE-117	8x3 in
AHWE-118	17x4 in

Urinal Stainless Steel (Cap. 700ml) Economy Model

AHW-121	Male
AHW-122	Female



Bed Pan Stainless Steel

AHWE-121	Perfection without cover
AHWE-122	Perfection with cover
AHWE-123	Male with cover

Irrigator (Douch Can) Stainless Steel Economy

AHW-128	With Handle, Hook & Nipple 2 Ltr.
AHW-129	With Handle, Hook & Nipple 1.5 Ltr.





Gallipot, Stainless Steel

- AHW-147 50ml., without Lid
- AHW-148 100ml., without Lid
- AHW-149 150ml., with Lid
- AHW-150 250ml., with Lid

Foreceps Jar Stainless Steel

- AHW-151 With Heavy Base, 2x7.25 in
- AHW-152 With Heavy Base, 2x7.25 in



Thermometer Jar, Stainless Steel

AHW-153

Dressing Drum, Stainless Steel

- AHW-103 225 X 225 mm (9 x 9 in.)
- AHW-104 275 X 240 mm (11 x 9.5 in.)
- AHW-105 275 X 125 mm (11 x 5 in.)
- AHW-106 350 X 240 mm (14 x 9.5 in.)
- AHW-107 350 X 130 mm (14 x 5.2 in.)
- AHWE-100 Size 9 x 9 in.
- AHWE-101 Size 11 x 9 in.
- AHWE-102 Size 14 x 9 in.



Miscellaneous Products



Medicine Cup



Iodine Cup



Dustbin



Measuring Jug



Badpan (Plastic)



Forceps Jar (Plastic)



Kidney Tray (Plastic)



Measuring Jug (Plastic)



Urinal Pot (Plastic)



Thermometer Jar (Plastic)

A TRUE PASSION FOR BETTER HEALTHCARE



WorldPrime
EXPORT LLP

MileStone, Near Zanzarda Cross Road, Junagadh 362001. Gujarat (INDIA)
Cell. +91 8866197888, +91 8866196888
sales@worldprimeexport.com