



“AN
ISO 9001- 2015
WHO - GMP
CERTIFIED CO.”



- * DIAGNOSTIC KITS**
- * LAB CHEMICALS**
- * LAB REAGNTS**
- * LAB GLASS WARES**
- * LAB ACCESSORIES**

CHEMISTRY LIQUID REAGENTS

PRODUCTS	METHODS	WORKING REAGENT RATIO	PACK SIZE(ml)
Albumin	BCG	Single Reagent	2x50
Albumin	BCG	Single Reagent	2x125
Alkaline Phosphatase	p-NPP with DEA Buffer	2 Reagent (4:1 Ratio)	25
Alkaline Phosphatase	p-NPP with DEA Buffer	2 Reagent (4:1 Ratio)	125
α-Amylase	GAL G2-CNP	Single Reagent	1x10
α-Amylase	GAL G2-CNP	Single Reagent	1x50
α-Amylase	GAL G2-CNP	Single Reagent	2x50
SGOT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	1x25
SGOT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	125
SGOT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	5x50
SGOT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	500
SGPT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	1x25
SGPT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	125
SGPT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	5x50
SGPT (UV Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	500
Bilirubin Total &Direct	Diazo - End Point	2 Reagent (50:1 Ratio)	2x50
Bilirubin Total &Direct	Diazo - End Point	2 Reagent (50:1 Ratio)	2x125
Bilirubin Total &Direct	Diazo - End Point	2 Reagent (50:1 Ratio)	500
Bilirubin Direct	Diazo - End Point	2 Reagent (50:1 Ratio)	2x50
Bilirubin Direct	Diazo - End Point	2 Reagent (50:1 Ratio)	2x125
Bilirubin Direct	Diazo - End Point	2 Reagent (50:1 Ratio)	500
Bilirubin Total	Diazo - End Point	2 Reagent (50:1 Ratio)	2x50
Bilirubin Total	Diazo - End Point	2 Reagent (50:1 Ratio)	2x125
Bilirubin Total	Diazo - End Point	2 Reagent (50:1 Ratio)	500
Calcium	O-cpc	2 Reagent (1:1 Ratio)	2x25
Calcium	O-cpc	2 Reagent (1:1 Ratio)	2x50
Calcium	O-cpc	2 Reagent (1:1 Ratio)	2x125
Cholesterol + HDL PPT	CHOD-POD, Phosphotungstate	Single Reagent	2x50
Cholesterol	CHOD - POD	Single Reagent	1x50
Cholesterol	CHOD - POD	Single Reagent	2x50
Cholesterol	CHOD - POD	Single Reagent	500
CK-Nac (Kinetic)	ACTIVATED NAC	2 Reagent (4:1 Ratio)	25
CK-Nac (Kinetic)	ACTIVATED NAC	2 Reagent (4:1 Ratio)	125
Chloride	Mercuric thiocyanate	Single Reagent	2x20
Chloride	Mercuric thiocyanate	Single Reagent	2x125
CK-MB (Kinetic)	Immunoinhibition	2 Reagent (4:1 Ratio)	25
CK-MB (Kinetic)	Immunoinhibition	2 Reagent (4:1 Ratio)	125
Cleaning Solution (Semi Auto Analyser)			2x125
Creatinine (Jaffe's 2P Kinetic)	Alkaline Picrate	2 Reagent (1:1 Ratio)	2x50
Creatinine (Jaffe's 2P Kinetic)	Alkaline Picrate	2 Reagent (1:1 Ratio)	2x125
Creatinine (Jaffe's 2P Kinetic)	Alkaline Picrate	2 Reagent (1:1 Ratio)	4x125
GGT (Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	25
GGT (Kinetic)	Optimized IFCC	2 Reagent (4:1 Ratio)	125
Glucose	GOD-POD	Single Reagent	1000
Glucose	GOD-POD	Single Reagent	2x125
LDH (Kinetic)	IFCC	2 Reagent (4:1 Ratio)	25
LDH (Kinetic)	IFCC	2 Reagent (4:1 Ratio)	125
Micro Protein	Pyrogallol red	Single Reagent	2X50
Micro Protein	Pyrogallol red	Single Reagent	2x125
Phosphorus	UV-End point	Single Reagent	2X50
Phosphorus	UV-End point	Single Reagent	2x125

CHEMISTRY LIQUID REAGENTS

PRODUCTS	METHODS	WORKING REAGENT RATIO	PACK SIZE(ml)
Triglycerides	GPO-POD	Single Reagent	1x50
Triglycerides	GPO-POD	Single Reagent	2x50
Triglycerides	GPO-POD	Single Reagent	10x50
Total Protein	Biuret Method	Single Reagent	2x50
Total Protein	Biuret Method	Single Reagent	2x125
Uric Acid	Uricase-POD	Single Reagent	1x50
Uric Acid	Uricase-POD	Single Reagent	2x50
Uric Acid	Uricase-POD	Single Reagent	10x50
Urea UV (2P Kinetic)	Urease-GLDH	2 Reagent (4:1 Ratio)	25
Urea UV (2P Kinetic)	Urease-GLDH	2 Reagent (4:1 Ratio)	125
Acid Phosphatase (kinetic)	PNPP	2 Reagent (4:1 Ratio)	25
Acid Phosphatase (kinetic)	PNPP	2 Reagent (4:1 Ratio)	125
Cholesterol Direct(EP), HDL	CHOD-POD, Phosphotungstate	2 Reagent (3:1 Ratio)	20
Cholesterol Direct(EP), HDL	CHOD-POD, Phosphotungstate	2 Reagent (3:1 Ratio)	80
Cholesterol Direct(EP), LDL	Direct Enzymatic Chlorimetric	2 Reagent (3:1 Ratio)	20
Cholesterol Direct(EP), LDL	Direct Enzymatic Chlorimetric	2 Reagent (3:1 Ratio)	80
Lipase (Kinetic)	Enzymatic	2 Reagent (5:1 Ratio)	1x12
Lipase (Kinetic)	Enzymatic	2 Reagent (5:1 Ratio)	4x12

CLINICAL CHEMISTRY (LYOPHILIZED/POWDER)

PRODUCTS	METHOD	WORKING REAGENT RATIO	PACK SIZE(ml)
SGPT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	4x20
SGPT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	2x50
SGPT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	5x50
SGPT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	20x50
SGOT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	4x20
SGOT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	2x50
SGOT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	5x50
SGOT (UV Kinetic 30 Days)	Optimized IFCC	Single Reagent	20x50
Cholesterol (EP 90 Days)	CHOD-POD	Single Reagent	2x50
Cholesterol (EP 90 Days)	CHOD-POD	Single Reagent	5x50
Cholesterol (EP 90 Days)	CHOD-POD	Single Reagent	20x50
Glucose (EP 60 Days)	GOD-POD	Single Reagent	4x125
Triglycerides (EP 30 Days)	GPO-POD	Single Reagent	5x10
Triglycerides (EP 30 Days)	GPO-POD	Single Reagent	2x50
Triglycerides (EP 30 Days)	GPO-POD	Single Reagent	5x50
Uric Acid (EP 30 Days)	Uricase-POD	Single Reagent	2x50
Uric Acid (EP 30 Days)	Uricase-POD	Single Reagent	5x50
Uric Acid (EP 30 Days)	Uricase-POD	Single Reagent	20x50
Urea Berthelot (EP 21 Days)	Berthelot	2 Reagent (1:1 Ratio)	2x50
Urea UV (2P Kinetic 30 Days)	Urease-GLDH	Single Reagent	4x20
Urea UV (2P Kinetic 30 Days)	Urease-GLDH	Single Reagent	5x50
Urea UV (2P Kinetic 30 Days)	Urease-GLDH	Single Reagent	20x50

SEROLOGY - AGGLUTINATION

PRODUCTS	METHOD	PACK SIZE
ASO	Latex Slide Agglutination	25T
ASO	Latex Slide Agglutination	50T
ASO	Latex Slide Agglutination	100T
CRP	Latex Slide Agglutination	25T
CRP	Latex Slide Agglutination	50T
CRP	Latex Slide Agglutination	100T
RF	Latex Slide Agglutination	25T
RF	Latex Slide Agglutination	50T
RF	Latex Slide Agglutination	100T
RPR	Carbon Slide Flocculation	50T
RPR	Carbon Slide Flocculation	100T
RPR	Carbon Slide Flocculation	500T
Widal	Slide Agglutination	4x5ml
Widal - OH	Slide Agglutination	(2+2)x5ml

RAPID TEST KITS

PRODUCTS	SAMPLE REQUIREMENT	TECHNOLOGY	PACK SIZE
Dengue IgG/IgM CARD	Serum/Plasma/W.B	Lateral Flow Chromatography	10T
Dengue Ns1 CARD	Serum/Plasma/W.B	Lateral Flow Chromatography	10T
Dengue Ns1/IgG/IgM Combo	Serum/Plasma/W.B	Lateral Flow Chromatography	10T
Malaria PF/PV Antigen	Whole Blood	Lateral Flow Chromatography	40T
Malaria PF/Pan Antigen	Whole Blood	Lateral Flow Chromatography	40T
Preg. Midstream Card- Clear Response Plus	Urine	Lateral Flow Chromatography	01T
Preg. Card- Clear Response	Urine	Lateral Flow Chromatography	40T
Preg. Card - Clear & Sure	Urine	Lateral Flow Chromatography	50T
Preg. Card - Clear Sign Plus	Urine	Lateral Flow Chromatography	50T
Preg. Strip - Clear Sign	Urine	Lateral Flow Chromatography	100T
FSH Test Card - Clear Menopase	Urine	Lateral Flow Chromatography	3T
LH Test Card - Clear Ovum LH	Urine	Lateral Flow Chromatography	10T
Fertility Test - Clear Ovum	Urine	Lateral Flow Chromatography	5+1T
Syphilis Card	Serum	Lateral Flow Chromatography	50T
Syphilis Strip	Serum	Lateral Flow Chromatography	100T
S.Typhi IgG/IgM Card	Serum/Plasma	Lateral Flow Chromatography	40T
Troponin- I Card	Serum/Plasma/W.B	Lateral Flow Chromatography	10T
UTI Test	Urine	Lateral Flow Chromatography	01T
Male Fertility Test Kit (Sperm Test)	Urine	Lateral Flow Chromatography	5T
Anti Ovarian Antibodies Test Kit	Urine	Lateral Flow Chromatography	01T

VACUTAINERS GLASS, 13X75

DESCRIPTION	PACK SIZE(Pcs)
PLAIN, Vacutainer Tubes, 4ml Red	100
EDTA K2(K3), 2ml Purple	100
EDTA NAF, 2ml, Grey	100
Sodium Heparin, 4ml, Green	100
Sodium Citrate 3.2%, 2ml Blue	100
Sodium Citrate 3.8%, 2ml Black	100
* Gel & Clot Activator, 5ml Yellow	100

MICROPIPETTE TIPS

DESCRIPTION	PACK SIZE(Pcs)
Blue Tips Gilson type 1000ul	500
White Tips Gilson type 1000ul	500
Yellow Tips, Gilson type 200ul	1000
White Tips Gilson type 200ul	1000
White Tips Gilson type 10 ul	1000
White Tips 1000ul with filter	500
White Tips 200ul with filter	1000
White Tips 10ul with filter	1000
White Tips 1000ul RNA/DNA Free	500
White Tips 200ul RNA/DNA Free	1000
White Tips 10u RNA/DNA Free	1000
White Tips 1000ul (with filter) in Rack	60
White Tips 200ul (with filter) in Rack	96
White Tips 10ul (with filter) in Rack	96

LAB GLASS WARES

DESCRIPTION	PACK SIZE
Cover Glass Square 18x18mm (Indian)	20x10gm
Cover Glass Square 22x22mm (Indian)	20x10gm
Microslides 75mmx25mmx1.25mm(Indian)	10x50Nos
Microslides 75mmx25mmx1.30mm(Indian)	10x50Nos
ESR Pipette Westergreen(Glass)	100
Haemometer Pipette 20µl(Glass)	100
Haemometer Tube Round(Glass)	100
HB Square Tube(Glass)	10
RBC Pipette Pvc Tubing(Glass)	100
WBC Pipette Pvc Tubing(Glass)	100
Winthrope Tube(Glass)	100
CoverGlasses, 18x18mm,100pcs, #1 (Imported)	10
CoverGlasses, 22x22mm 100pcs, #1 (Imported)	10
Cover glass 22 x 50 mm 100Pcs, # 1(Imported)	5
Microslides 75mmx25mmx1.25mm(Imported)	50
Microslides 75mmx25mmx1.25mm(Imported) Frosted one end	50
Test Tube Glass Disposable 12x75mm	100
Test Tube Glass Disposable 12x100mm	100
Haematocrit Capillaries Heparinised	20x100
Haematocrit Capillaries Non Heparinised	20x100

LAB PLASTIC WARES

DESCRIPTION	PACK SIZE(Pcs)
* Storage Vial 2ml	100
Petridish 90mm ETO Sterilized	100
Urine Container 30ml ETO-Sterilized	100
Urine Container 50ml ETO-Sterilized	100
Urine Container 60ml ETO-Sterilized	100
Urine Container 120 ml	50
Stool Container 30 ml PP	50
Stool Container 30 ml PS	50
Centrifuge Tube 15ml(Sterilized)	100
Centrifuge Tube 50ml(Sterilized)	50
Pasteur Pipettes 3ml(Sterilized)	500
Test Tube (Ria Tube) 12x75mm PS	500
Microcentrifuge Tube 0.5ml conical (DNA/RNA FREE)	1000
Microcentrifuge Tube 1.5ml conical (DNA/RNA FREE)	500
Microcentrifuge Tube 2.0ml conical (DNA/RNA FREE)	500
Microcentrifuge Tube 0.5ml conical	1000
Microcentrifuge Tube 1.5ml conical	500
Microcentrifuge Tube 2.0ml conical	500
Centrifuge Tube 15ml Conical	100
Centrifuge Tube 50 ml conical	50
Cryovial Tube 1.8ml	500
Sample Cups Hitachi	500
Pasteur Pipettes 3ml	500
* Serological Pipette 1ml	500
* Serological Pipette 2ml	500
* Serological Pipette 5ml	200
* Serological Pipette 10ml	200
Serological Pipette 25ml	150

MICROPIPETTES

DESCRIPTION	PACK SIZE(Pcs)
Micropipette Fixed 10 μ L - Indian	1
Micropipette Fixed 20 μ L - Indian	1
Micropipette Fixed 100 μ L - Indian	1
Micropipette Fixed 500 μ L - Indian	1
Micropipette Fixed 1000 μ L - Indian	1
Micropipette Variable 5-50 μ L - Indian	1
Micropipette Variable 10-100 μ L - Indian	1
Micropipette Variable 100-1000 μ L - Indian	1
Micropipette holder stand - Indian	1
Micropipette Variable 5-50 μ L - Imported	1
Micropipette Variable 10-100 μ L - Imported	1
Micropipette Variable 100-1000 μ L - Imported	1

TOURNIQUETES

DESCRIPTION	PACK SIZE(Pcs)
Tourniquete Buckle Premium	100
Tourinquete Standard	100

LAB ACCESSORIES

DESCRIPTION	PACK SIZE(Pcs)
Test Tube Rack D13 mm, 50 Holes	10
Test Tube Rack D18 mm, 50 Holes	10
Test Tube Rack 50 wells Multipurpose	10
Rack for Gilson Tips 200µl 96 Hole	10
Rack for Gilson Tips 1000µl 60 Hole	10
Counting Chambers	10
ESR Pipette Disposable	200
Blood Cell Counter 5 Keys	1
Blood Cell Counter 8 Keys	1
Alcohol Swabs Saturated	100
Embedding Cassettes	100
Embedding Ring	100
Pipette Pump 2ml	1

SLIDE BOXES

DESCRIPTION	PACK SIZE(Pcs)
Slide Box for 1 Pcs	100
Slide Box for 2 Pcs	100
Slide Mailer for 5 Slides	100
Slide Box for 25 Slides	1
Slide Box for 50 Slides	1
Slide Box for 100 Slides	17

THERMAL PAPERS

DESCRIPTION	PACK SIZE(Pcs)
Thermal Paper 110mmx25mtrs	1
Thermal Paper 57mmx25mtrs	1
Thermal Paper 80mmx25mtrs	1
Thermal Paper 50mmx25mtrs	1

BLOOD LANCETS

DESCRIPTION	PACK SIZE(Pcs)
Blood Lancets, Steel Type	200
Blood Lancets Plastic Screw	100
Blood Lancets Plastic Flat	100

GYNAECOLOGICAL RANGE

DESCRIPTION	PACK SIZE(Pcs)
Cytobrush, standard type, nylon bristle,ss wire, PP stick	100
Swab Sticks Wooden 6 inch	100
Swab Sticks Sterilized Wooden (in plastic tubes)	100
Swab Sticks Sterilized Plastic(in plastic tubes)	100
Wooden Ayre Scrapper / Spatula	100
Papsmear Test Kit - 4pcs - Sterilized	100
Papsmear Test Kit - 2pcs - Sterilized	100
Swab Stick Sterilized	100
Endometrial Suction Curette	100
* Umbilical Cord Clamp - Paper Pouch	100
* Mucus Extractor - Plastic Pouch	1
Cotton Applicator	100

OTC PRODUCTS

PRODUCTS	PACK SIZE(Pcs)
Hot Water Bottle Electric Standard	1
Hot Water Bottle Electric Premium	1
Hot Water Bottle 2000ml without cover standard	1
Hot Water Bottle 2500ml without cover - Premium Quality	1
Analogue Weighing Scale	1
Glass Top Digital Weighing Scale - 6mm Standard	1
Glass Top Digital Weighing Scale - 6mm Premium	1
Body Fat analyzer	1
Aneroid Sphygmomanometer with Dual Head Stethoscope	1
Stethoscope	1
Digital Thermometers Fixed Type	1
Compressor nebulizer	1
Dust Proof Face Mask(Equivalent to N95)	1
Surgical Tape 1/2" 9.1 meter	24
Surgical Tape 1" 9.1 meter	12
Surgical Tape 2" 9.1 meter	6
Surgical Tape 1/2" 5.0 meter	24
Surgical Tape 1" 5.0 meter	12
Surgical Tape 2" 5.0 meter	6

LAB CHEMICALS LIQUID REAGENTS

PRODUCTS	PACK SIZE(ml)
Sodium Hypochlorite Solution 5%	5000ml
Acetone	500ml
Acetone	5000ml
Formaldehyde Solution, 37-40%	500ml
Formaldehyde Solution, 37-40%	5000ml
Hcl N/10	500ml
Hcl N/10	5000ml
Normal Saline	500ml
Deionized water	5000ml
Ceederwood Oil	30ml
Ceederwood Oil	500ml
Absolute Alcohol	500ml
Isopropyl Alcohol	500ml
Isopropyl Alcohol	5000ml
Solution for Lishman Stain(Methanol)	500ml
Solution for Lishman Stain(Methanol)	2500ml
Acetic Acid Glacial	500ml
Barium Chloride 10%	500ml

LAB CHEMICALS LIQUID REAGENTS

PRODUCTS	PACK SIZE(ml)
Benedict's Reagent	500ml
Benedict's Reagent	5000ml
Carbol Fuchsin's Strong	500ml
CSF D.F	500ml
Crystal Violet Sol Aqs.	500ml
Crystal Violet Sol Alc.	500ml
Drabkin's Sol. With STD	1000ml
Drabkin's Sol. With STD	5000ml
Ehrlich's Reagent	125ml
EDTA Disodium	100gm
EDTA Sol. 5%	500ml
EDTA N/10	500ml
Esbach's Rea.	500ml
Field Stain A	500ml
Field Stain B	500ml
Gentian Violet Aqueous Sol.	500ml
Giemsa's Sol.	500ml
Glycerine(Glycerol)	500ml
Leishman's Stain Sol.	500ml
Lysol	500ml
Lugol's Iodine	500ml
Laboclean Liq.Ditergent	5000ml
Methylene Blue Sol.Alk.	500ml
R.B.C. D.F.	500ml
Safranin Stain Sol.	500ml
Semen D.F.	500ml
Sodium Chloride Sol. 0.85%	500ml
Sodium Citrate 3.8%	500ml
W.B.C. D.F.	500ml
Xylene Pure	500ml
Immersion Oil	30ml
Immersion Oil	500ml

LAB CHEMICALS & REAGENT (POWDER)

PRODUCTS	PACK SIZE(gm)
* Crystal Violet Pow.	25gm
* Crystal Violet Pow.	100gm
* Giemsa's Pow.	25gm
* Giemsa's Pow.	100gm
* Leishman's Stain Pow.	25gm
* Leishman's Stain Pow.	100gm
* Methylene Blue	25gm
* Methylene Blue	100gm

OTHER PRODUCTS

PRODUCTS	PACK SIZE(Pcs)
* Waste Carry Bag - Bio Bag	50
Goggles - Plastic Transparent	10
* Blood Roller Mixer	1
Eyeshied	10

ID BRACELETS

PRODUCTS	PACK SIZE (Pcs)
ID bracelet for children,written type	100
ID bracelet for Adult,written type	100
ID bracelet for children,Card Insert type	100
ID bracelet for Adult,Card Insert type	100
ID bracelet for Adult & Child,Card Insert type	100

URINE STRIPS

PRODUCTS	PARAMETERS	PACK SIZE
Urisign 2P	Glucose, Protein	100 Strip Bottle
Urisign 3P	Glucose, Protein, Ketone	100 Strip Bottle
Urisign 10P	Nitrate, pH, SG, Blood, Glu, Bil, Urobilinogen, Ketone, leucocytes, Protein	100 Strip Bottle
Urisign 4P	Glucose, Ketone, Protein, pH	100 Strip Bottle

COVID-19 PRODUCT

PRODUCTS

PACK SIZE

VTM Kit	50Pcs
Nasal Swab	100 Pair
Oral Swab	100 Pair

SYRINGES

PRODUCTS

PACK SIZE(Pcs)

Syringe 1 ml	100 pcs
Syringe 2 ml	100 pcs
Syringe 5 ml	100 pcs
Syringe 10 ml	100 pcs
Syringe 20 ml	100 pcs



A TRUE PASSION FOR BETTER HEALTHCARE



WorldPrime
EXPORT LLP

MileStone, Near Zanzarda Cross Road, Junagadh 362001. Gujarat (INDIA)
Cell. +91 8866197888, +91 8866196888
Email:-worldprimeexport@gmail.com